# **PRODUCTINFO**

## INSTRUMENTS, CHEMICALS AND CONSUMABLES FOR CELL DIAGNOSTICS

## PRODUCT GROUP CHEMICALS

### SPRAYFIX®

Sprayfix® is applied immediately after smearing the cell specimen and before airdrying. Cells are protected with a tough, water-soluble film that preserves nuclear details and cellular morphology for microscopic examination. The smears are prepared in the conventional way. The slides dry in a few minutes and are ready to send to the laboratory in any standard mailer.

Sprayfix® is water-soluble and staining procedures may be commenced without first soaking in solvents. The 200-ml spray bottle is sufficient for processing apr. 1.500 slides.



200 ML ART NO: 01500N-ENG

#### **SPECIFICATIONS**

Shelflife: 3 years from manufacturing date in a sealed packaging..

Volume: 200 mL

Art code: Sprayfix 01500N-ENG



